

Recombinant Human IL7RA/CD127 Protein (His Tag)

RPES4708

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES4708

Purity: > 92 % as determined by reducing SDS-PAGE.

Contents: 100µg, 1mg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: CDW, IL7R, CD127, CDW127, IL-7R-alpha, IL7RA, ILRA, IL-7R, IL-7RA, IL-7R \hat{I} \pm , Interleukin-7 receptor subunit alpha, IL-7R&, alpha, CD 127, CD127 antigen, IL 7R, IL 7R alpha, IL-7 receptor subunit alpha, IL-7R subunit alpha, IL7Ralpha, Interleukin 7 receptor, Interleukin 7 receptor alpha chain, Interleukin 7 receptor isoform H5 6

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Measured by its ability to bind biotinylated mouse TSLP-his in functional ELISA.

Calculated MW: 26.3 kDa

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 45-50 kDa

Reconstitution: Please refer to the printed manual for detailed information.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Accession: NP_002176.2

Source: HEK293 Cells-derived Human IL7RA/CD127 protein Met 1-Gly 236, with an C-terminal His

Sequence: Met 1-Gly 236

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.