

Recombinant Human CD3 ϵ /CD3E Protein (Fc Tag)

RPES5124

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES5124

Purity: > 95 % as determined by reducing SDS-PAGE.

Contents: 50 μ g, 10 μ g
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: < 1.0 EU per μ g of the protein as determined by the LAL method.

Synonyms: CD3 epsilon, CD3E, CD3e, IMD18, T-Cell Surface Antigen T3/Leu-4 Epsilon Chain, T-Cell Surface Glycoprotein CD3 Epsilon Chain, T3E

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-Fc

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Immobilized Human Human CD3E-Fc at 2 μ g/ml (100 μ l/well) can bind Anti- Human/Monkey CD3E mAb. The ED50 of Anti-Human/Monkey CD3E mAb is 0.51 μ g/ml.

Calculated MW: 38.7 kDa

Formulation: Lyophilized from a 0.2 μ m filtered solution of PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

Observed MW: 40-50 kDa

Reconstitution: Please refer to the printed manual for detailed information.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Accession: P07766

Source: HEK293 Cells-derived Human CD3 epsilon, CD3E protein Asp23-Asp126, with an C- terminal Fc

Sequence: Asp23-Asp126

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.