

Recombinant Human IL-1 α Protein

RPES5954

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES5954

Purity: > 90% as determined by reducing SDS-PAGE.

Contents: 100 μ g, 20 μ g, 500 μ g, 1mg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: < 1.0 EU/mg of the protein as determined by the LAL method

Synonyms: IL1F, IL1A, IL-1A, IL1, IL1-ALPHA, IL1F1, BAF, FAF, Hematopoietin 1, Hematopoietin-1, IL 1 alpha, IL 1A, IL1 ALPHA, IL-1 alpha, IL1-A, Interleukin 1 alpha, Interleukin1 alpha, Interleukin-1 alpha, LAF, LEM, Preinterleukin 1 alpha, Pro interleukin 1 alpha, Pro-Interleukin-1-Alpha

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: HEK293 Cells

Bio-Activity: Not validated for activity

Calculated MW: 30.6 kDa

Formulation: Lyophilized from a 0.2 μ m filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Observed MW: 35 kDa

Reconstitution: It is recommended that sterile water be added to the vial to prepare a stock solution of 0.5 mg/mL. Concentration is measured by UV-Vis.

Accession: P01583

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Source: HEK293 Cells-derived Human IL-1 α protein Met1-Ala271, with an C-terminal His

Sequence: Met1-Ala271

Form: Lyophilized powder

this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.