

Recombinant Human BRD4 (N-10His-Flag)

RPES6330

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES6330

Purity: > 90 % as determined by reducing SDS-PAGE.

Contents: 50µg, 10µg
Bradford Reagent: 1 vial (2ml)

Concentration: Subject to label value.

Species: Human

Endotoxin: < 1.0 EU per µg of the protein as determined by the LAL method.

Synonyms: HUNK1, HUNK1bromodomain-containing protein 4, HUNKI, MCAP

Storage: Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles. Store Bradford Reagent at Room Temperature for 1 year.

Tag: N-His-Flag

Shipping: This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

Expression Host: E.coli

Bio-Activity: Not validated for activity

Calculated MW: 49.0 kDa

Formulation: Supplied as a 0.2 µm filtered solution of 50mM HEPES, 200mM NaCl, 1mM DTT, 10% Glycerol, pH 7.5.

Observed MW: 55-60 kDa

Reconstitution: -

Accession: O60885

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Source: E.coli-derived Human BRD4 protein Glu49-Glu460, with an N-terminal His & Flag

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Sequence: Glu49-Glu460

Notes: Centrifuge before opening to ensure complete recovery of vial contents.

Form: Liquid