

# Recombinant SARS-CoV-2 Spike Protein, Biotinylated (RBD, Avi-His RPES6808

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES6808

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 20µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** SARS-CoV-2

**Endotoxin:** < 1.0 EU per µg of the protein as determined by the LAL method.

**Synonyms:** NCP-CoV RBD, NCP-CoV Spike, NCP-CoV s1, NCP-CoV s2, RBD, S1, S2, Spike RBD, coronavirus spike, cov spike, ncov RBD, ncov s1, ncov s2, ncov spike, novel coronavirus RBD, novel coronavirus s1, novel coronavirus s2, novel coronavirus spike

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-Avi-His

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** HEK293 Cells

**Bio-Activity:** Recombinant Human ACE2 Protein (mFc Tag)(Active)(Cat#RPES0051) at 2 µg/mL (100 µL/well) can bind Recombinant 2019-nCoV Spike Protein, Biotinylated (RBD, Avi-His Tag)(Active)(Cat#RPES6808), the EC50 of RPES6808 is 40-80 ng/mL.

**Calculated MW:** 28.4 kDa

**Formulation:** Lyophilized from sterile PBS,pH 7.4 Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Observed MW:** -

**Reconstitution:** Please refer to the printed manual for detailed information.

**Accession:** YP\_009724390.1

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** HEK293 Cells-derived SARS-CoV-2 SARS-CoV-2 Spike(RBD) protein Arg319-Phe541, with an C-terminal Avi & His

**Sequence:** Arg319-Phe541

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Lyophilized powder