

# Recombinant SARS-CoV-2 Spike S1( $\Delta$ HV69-70, $\Delta$ Y144, N501Y, A570D, RPES7073

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES7073

**Purity:** > 90 % as determined by reducing SDS-PAGE.

**Contents:** 100 $\mu$ g  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** SARS-CoV-2

**Endotoxin:** < 1.0 EU per  $\mu$ g of the protein as determined by the LAL method.

**Synonyms:** 2019-nCoV S1 protein, S1 protein, coronavirus S1 Protein, cov S1 Protein

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** HEK293 Cells

**Bio-Activity:** Immobilized human ACE2 protein (Fc tag)(Cat: RPES1994) at 2  $\mu$ g/mL (100  $\mu$ L/well) can bind SARS-CoV-2 (2019-nCoV) Spike S1(HV69-70 deletion, Y144 deletion, N501Y, A570D, D614G, P681H)-His (RPES7073), the EC50 of RPES7073 is 10-60 ng/mL.

**Calculated MW:** 76.1 kDa

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Observed MW:** 112.2 kDa

**Reconstitution:** Please refer to the printed manual for detailed information. Immobilized human ACE2

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

protein (Fc tag)(Cat: RPES1994) at 2 µg/mL (100 µL/well) can bind SARS- CoV-2 (2019-nCoV) Spike S1(HV69-70 deletion, Y144 deletion, N501Y, A570D, D614G, P681H)-His (RPES7073), the EC50 of RPES7073 is 10-60 ng/mL.

**Accession:** YP\_009724390.1

**Source:** HEK293 Cells-derived SARS-CoV-2 SARS-CoV-2 Spike S1( $\Delta$ HV69-70,  $\Delta$ Y144, N501Y, A570D, D614G, P681H) protein Met1-Arg685(HV69-70 deletion, Y144 deletion, N501Y, A570D, D614G, P681H), with an C-terminal His

**Sequence:** Met1-Arg685(HV69-70 deletion, Y144 deletion, N501Y, A570D, D614G, P681H)

**Form:** Lyophilized powder

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

