

# Recombinant SARS-CoV-2 N Protein (S194L)(His Tag)

RPES7081

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES7081

**Purity:** > 85 % as determined by reducing SDS-PAGE.

**Contents:** 100µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** SARS-CoV-2

**Endotoxin:** Please contact us for more information.

**Synonyms:** SARS-CoV-2 coronavirus NP Protein, SARS-CoV-2 novel coronavirus, SARS-CoV-2 np Protein

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-His

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 46.6 kDa

**Formulation:** Lyophilized from sterile 50 mM PB, 500 mM NaCl, pH 7.0. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.

**Observed MW:** -

**Reconstitution:** Please refer to the printed manual for detailed information.

**Accession:** YP\_009724397.2

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in

**Source:** E.coli-derived SARS-CoV-2 SARS-CoV-2 N protein Met1-Ala419, with an N-terminal His

**Sequence:** Met1-Ala419

**Form:** Lyophilized powder

this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.