

# Recombinant SARS-CoV-2 Spike RBD (R408I)(His Tag)

RPES7131

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES7131

**Purity:** > 95 % as determined by reducing SDS-PAGE.

**Contents:** 500µg, 50µg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** Subject to label value.

**Species:** SARS-CoV-2

**Endotoxin:** Please contact us for more information.

**Synonyms:** 2019-nCoV Spike RBD Protein,  
2019-nCov RBD Protein

**Storage:** Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** C-His

**Shipping:** This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.

**Expression Host:** HEK293 Cells

**Bio-Activity:** Antigen was assessed by binding ability with ACE-2(RPES0028) or Anti-2019- nCoV-S1 mAb(5D9).

**Calculated MW:** 25.9 kDa

**Formulation:** Supplied as a 0.2 µm filtered solution of PBS, pH7.4.

**Observed MW:** 35 kDa

**Reconstitution:** -

**Accession:** QHD43416.1

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** HEK293 Cells-derived SARS-CoV-2 SARS-CoV-2 Spike RBD (R408I) protein Arg319- Phe541(R408I), with an C-terminal His

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Sequence:** Arg319-Phe541(R408I)

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Liquid