

RPES8211

Product Information

Product SKU:	RPES8211	Expression Host:	E.coli	Size:	20µg
Tag:	N-Sumo	Reactivity:	Mouse	Accession:	O88935-2

Additional Information

Calculated MW:	29.1 kDa	Observed MW:	37 kDa
Sequence:	Ser510-Pro656		

Protein Information

Background: -

Synonyms: -

Endotoxin: < 10 EU/mg of the protein as determined by the LAL method

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Purity: > 90% as determined by reducing SDS-PAGE.

Bio-Activity: Not validated for activity

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.