

Recombinant Pennsylvania firefly Luciferase Protein (His Tag)

RPES8217



Product Information

Product SKU:	RPES8217	Expression Host:	E.coli	Size:	20µg
Tag:	N-His	Reactivity:	Pennsylvania firefly	Accession:	Q27757

Additional Information

Calculated MW:	59.8 kDa	Observed MW:	60 kDa
Sequence:	Met1-Leu545		

Protein Information

Background: -

Synonyms: -

Endotoxin: < 10 EU/mg of the protein as determined by the LAL method

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Purity: > 90% as determined by reducing SDS-PAGE.

Bio-Activity: Not validated for activity

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.