

# Recombinant Human S100A8 Protein(GST Tag)

RPES9358

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES9358

**Purity:** > 95% as determined by reducing SDS-PAGE.

**Contents:** 100µg, 20µg, 500µg, 1mg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Human

**Endotoxin:** < 10 EU/mg of the protein as determined by the LAL method

**Synonyms:** Protein S100-A, S100a, Kcnip, KCHIP, MRP, S100A8, 60B8AG, CAGA, CFAG, CGLA, CP-10, L1Ag, MA387, MIF, MRP8, NIF, P8, Calgranulin A, Calgranulin-A, Cystic fibrosis antigen, Leukocyte L1 complex light chain, MRP-8, Protein S100-A8, CALP, KCNIP4, KCHIP4, AI323541, B8Ag, BEE11, Calprotectin, Calprotectin L1L subunit, Chemotactic cytokine CP-10, included, Migration inhibitory factor-related protein 8, Myeloid-related protein 8, Neutrophil cytosolic 7 kDa protein, Pro-inflammatory S100 cytokine, S100 calcium binding protein A8, S100 calcium binding protein A8 (calgranulin A), S100 calcium-binding protein A8, S100A8/S100A9 complex, S10A8, Urinary stone protein band A

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-GST

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 36.2 kDa

**Formulation:** Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Observed MW:** 40 kDa

**Reconstitution:** It is recommended that sterile water be added to the vial to prepare a stock solution of 0.5 mg/mL. Concentration is measured by UV-Vis.

**Accession:** P05109

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived Human S100A8 protein Met1-Glu93, with an N-terminal GST

**Sequence:** Met1-Glu93

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Lyophilized powder