

# Recombinant Mouse IL-33 Protein(Trx Tag)

RPES9380

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES9380

**Purity:** > 95% as determined by reducing SDS-PAGE.

**Contents:** 100µg, 20µg, 500µg, 1mg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Mouse

**Endotoxin:** < 10 EU/mg of the protein as determined by the LAL method

**Synonyms:** DVS, NF-HEV and IL- 1F, Interleukin-1 family member, C9orf, C9orf26, DVS27, Interleukin-1 family member 11, NF-HEV and IL- 1F11, NKHEV, IL33, Interleukin-33, IL-33, Interleukin 33, IL1F11, NFEHEV, NF-HEV

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.  
Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-Trx

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 37.4 kDa

**Formulation:** Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

**Observed MW:** 37 kDa

**Reconstitution:** It is recommended that sterile water be added to the vial to prepare a stock solution of 0.5 mg/mL. Concentration is measured by UV-Vis.

**Accession:** Q8BVZ5

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Source:** E.coli-derived Mouse IL-33  
protein Gln109-ILE266, with an N-terminal  
Trx

**Sequence:** Gln109-Ile266

**Form:** Lyophilized powder

**Notes:** Centrifuge before opening to ensure complete  
recovery of vial contents.