

# Recombinant Human S100A9 Protein(Trx Tag)

RPES9401

## Description

---

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Protein Information

---

**SKU:** RPES9401

**Purity:** > 95% as determined by reducing SDS-PAGE.

**Contents:** 100µg, 20µg, 500µg, 1mg  
Bradford Reagent: 1 vial (2ml)

**Concentration:** -

**Species:** Human

**Endotoxin:** < 10 EU/mg of the protein as determined by the LAL method

**Synonyms:** MAC, Protein S100-A, S100 calcium-binding protein A, Migration inhibitory factor-related protein, S100A, MRP, S100A9, 60B8AG, CAGB, CFAG, CGLB, L1AG, LIAG, MAC387, MIF, MRP14, NIF, P14, Calgranulin-B, Calprotectin L1H subunit, Leukocyte L1 complex heavy chain, Migration inhibitory factor-related protein 14, MRP-14, Protein S100-A9, S100 calcium-binding protein A9, Calgranulin B, Cystic fibrosis antigen B, Migration inhibitory factor related protein 14, MRP 14, Myeloid-related protein 14, OTTHUMP00000015331, S100 A9, S100 calcium binding protein A9, S100 calcium binding protein A9 calgranulin B, S10A9

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

**Tag:** N-Trx

**Shipping:** This product is provided as lyophilized powder which is shipped with ice packs.

**Expression Host:** E.coli

**Bio-Activity:** Not validated for activity

**Calculated MW:** 32.5 kDa

**Formulation:** Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

**Manufacturers Statement:** This final kit system is assembled and quality-released by Assay Genie Limited.

**Observed MW:** 30 kDa

**Reconstitution:** It is recommended that sterile water be added to the vial to prepare a stock solution of 0.5 mg/mL. Concentration is measured by UV-Vis.

**Accession:** P06702

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

**Source:** E.coli-derived Human S100A9 protein Met1-Pro114, with an N-terminal Trx

**Sequence:** Met1-Pro114

**Notes:** Centrifuge before opening to ensure complete recovery of vial contents.

**Form:** Lyophilized powder