

Recombinant Human IL17RA Protein(His Tag)

RPES9416

Description

This high-purity recombinant protein is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Protein Information

SKU: RPES9416

Purity: > 95% as determined by reducing SDS-PAGE.

Contents: 100µg, 20µg, 500µg, 1mg
Bradford Reagent: 1 vial (2ml)

Concentration: -

Species: Human

Endotoxin: < 1.0 EU/mg of the protein as determined by the LAL method

Synonyms: IMD, MGC, CTLA, CANDF, CDw, IL17RA, CANDF5, CD217, CDw217, IL-17RA, IL17R, IMD51, hIL-17R, CD127, CD217 antigen, CDW127, CTLA8, HIL 17R, hIL17R, I17RA, IL 17 receptor, IL 17 receptor A, IL 17RA, IL17, IL-17 R alpha, IL-17 receptor A, IL17A receptor, IL-17R, IL17R-A, IL17-RA, IL-17R-alpha, IL7RA, ILRA, Interleukin 17 receptor A, Interleukin-17 receptor A, MGC10262, CDw217, Interleukin-17 receptor A

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80 °C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Store Bradford Reagent at Room Temperature for 1 year.

Tag: C-His

Shipping: This product is provided as lyophilized powder which is shipped with ice packs.

Expression Host: Mammalian

Bio-Activity: Not validated for activity

Calculated MW: 31.6 kDa

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Observed MW: 50-65 kDa

Reconstitution: It is recommended that sterile water be added to the vial to prepare a stock solution of 0.5 mg/mL. Concentration is measured by UV-Vis.

Manufacturers Statement: This final kit system is assembled and quality-released by Assay Genie Limited.

Accession: Q96F46

Source: Mammalian-derived Human IL17RA protein Leu33-Trp320, with an C-terminal His

Sequence: Leu33-Trp320

Form: Lyophilized powder

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol

Notes: Centrifuge before opening to ensure complete recovery of vial contents.