

HDBS0166

Product Information

Product SKU:	HDBS0166	Host Species:	Humanized	Size:	100 µg
Target:	AXL	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	P30530
Isotype:	IgG1(E239D, M241L)		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	UFO
Common Name:	BA301, F1-2-15
Recommended Dilutions:	ELISA 1:5000-10000, Flow Cyt 1:100
Applications:	ELISA, Flow Cyt
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.