

# Anti-CD79B (Polatuzumab biosimilar) antibody, Rabbit IgG Chimeric mAb HDBS0207



## Product Information

<b>Product SKU:</b>	HDBS0207	<b>Host Species:</b>	Chimeric Rabbit	<b>Size:</b>	100 µg
<b>Target:</b>	CD79B	<b>Reactivity:</b>	Human		

## Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot:</b>	P40259
<b>Isotype:</b>	IgG		

## Protein Information

<b>Background:</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Synonyms:</b>	AGM6;B29;IGB
<b>Common Name:</b>	RO5541077、RG-7596
<b>Recommended Dilutions:</b>	ELISA 1:5000-10000, Flow Cyt 1:100
<b>Applications:</b>	ELISA, Flow Cyt
<b>Usage:</b>	Research use only
<b>Formulation:</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage:</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.