

**Biotinylated Anti-TIGIT (Tiragolumab biosimilar)
mAb**

HDBS0242



Product Information

Product SKU:	HDBS0242	Host Species:	Homo sapiens	Size:	100 µg
Target:	TIGIT	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	Q495A1
Isotype:	IgG1		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	TIGIT;SIG9;STM3
Common Name:	MTIG-7192-A, MTIG-7192A, MTIG7192A, RG-6058, RO-7092284, RO7092284
Recommended Dilutions:	ELISA 1:5000-10000; Flow Cyt 1:100
Applications:	ELISA, Flow Cyt
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.