

**Biotinylated Anti-OX40 (ivuxolimab
biosimilar)mAb**
HDBS0263



Product Information

Product SKU:	HDBS0263	Host Species:	Homo sapiens	Size:	100 µg
Target:	OX40	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P43489
Isotype:	IgG2		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	TNFRSF4;OX40;CD134;OX40L receptor;ACT35;TXGP1L
Common Name:	PF-04518600
Recommended Dilutions:	Flow Cyt 1:100
Applications:	Flow Cyt
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.