

**Biotinylated Anti-TNFSF11 (Denosumab
biosimilar)mAb**
HDBS0271



Product Information

| | | | | | |
|---------------------|----------|----------------------|--------------|--------------|--------|
| Product SKU: | HDBS0271 | Host Species: | Homo sapiens | Size: | 100 µg |
| Target: | TNFSF11 | Reactivity: | Human | | |

Additional Information

| | | | |
|-------------------|--------------|-----------------|--------|
| Conjugate: | Biotinylated | Uniprot: | O14788 |
| Isotype: | IgG2 | | |

Protein Information

| | |
|-------------------------------|---|
| Background: | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Synonyms: | CD254;hRANKL2;ODF;OPGL;OPTB2;RANKL;sOdf;TNLG6B;TRANCE |
| Common Name: | AMG162 |
| Recommended Dilutions: | Flow Cyt 1:100 |
| Applications: | Flow Cyt |
| Usage: | Research use only |
| Formulation: | Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Storage: | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |