

Biotinylated Anti-Her2 (Trastuzumab biosimilar) mAb

HDBS0277



Product Information

Product SKU:	HDBS0277	Host Species:	Humanized	Size:	100 µg
Target:	Her2	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P04626
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: ERBB2;CD340;HER-2/neu;HER2;MLN19;NEU;NGL;TKR1

Common Name: 4D5-8,4D5V8,Herceptin,rhuMabHER2

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.