

Anti-CCR2 (Plozalizumab Biosimilar) Biotinylated Antibody

HDBS0293

Description

This Anti-CCR2 (Plozalizumab Biosimilar) Biotinylated Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDBS0293
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Category:	Biosimilar Antibody
Synonyms:	CC-CKR-2, CCR2A, CCR2B, CD192, CKR2, CKR2A, CKR2B, CMKBR2, MCP-1-R
Applications:	Flow Cyt
Conjugation:	Biotinylated
Reactivity:	Human
Common Name:	MLN-1202, MLN1202, anti-CCR2, hu1D9

Antibody Data

Isotype:	IgG1
Uniprot ID:	P41597
Host Species:	Humanized
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Target:	CCR2

Preparation & Storage

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions: Flow Cyt 1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol