

Biotinylated Anti-KLRG1 (Ulviprubart biosimilar)
mAb

HDBS0355



Product Information

Product SKU:	HDBS0355	Host Species:	Humanized	Size:	100 µg
Target:	KLRG1	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	Q96E93
Isotype:	IgG1(E239D, M241L)		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	2F1;CLEC15A;MAFA;MAFA-2F1;MAFA-L;MAFA-LIKE
Common Name:	ABC-008,ABC008
Recommended Dilutions:	ELISA 1:5000-10000, Flow Cyt 1:100
Applications:	ELISA, Flow Cyt
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.