

Biotinylated Anti-FOLR1 (Farletuzumab biosimilar) mAb

HDBS0356



Product Information

Product SKU:	HDBS0356	Host Species:	Humanized	Size:	100 µg
Target:	FOLR1	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P15328
Isotype:	IgG1(E239D, M241L)		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: FBP;FOLR;FRalpha;NCFTD

Common Name: M3,MORAb-003

Recommended Dilutions: ELISA 1:5000-10000, Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.