

**Biotinylated Anti-CEACAM5 (Tusamitamab
biosimilar) mAb**
HDBS0388



Product Information

Product SKU:	HDBS0388	Host Species:	Chimeric	Size:	100 µg
Target:	CEACAM5	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P06731
Isotype:	IgG1(E239D, M241L)		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	CEA
Common Name:	SAR-408377
Recommended Dilutions:	ELISA 1:5000-10000, Flow Cyt 1:100
Applications:	ELISA, Flow Cyt
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.