

**Biotinylated Anti-GM-CSFR (Mavrilimumab biosimilar) mAb**  
**HDBS0408**



**Product Information**

<b>Product SKU:</b>	HDBS0408	<b>Host Species:</b>	Homo sapiens	<b>Size:</b>	100 µg
<b>Target:</b>	GM-CSFR	<b>Reactivity:</b>	Human		

**Additional Information**

<b>Conjugate:</b>	Biotinylated	<b>Uniprot:</b>	P15509
<b>Isotype:</b>	IgG4(P108S)		

**Protein Information**

<b>Background:</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Synonyms:</b>	GM-CSF-R-alpha; GMCSFR-alpha; GMR-alpha; CDw116; CD116; CSF2RA; CSF2R; CSF2RY
<b>Common Name:</b>	CAM-3001
<b>Recommended Dilutions:</b>	Flow Cyt 1:100
<b>Applications:</b>	Flow Cyt
<b>Usage:</b>	Research use only
<b>Formulation:</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage:</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.