

**Biotinylated Anti-IL12B (Ebdarokimab biosimilar)  
mAb**

**HDBS0418**



**Product Information**

<b>Product SKU:</b>	HDBS0418	<b>Host Species:</b>	Humanized	<b>Size:</b>	100 µg
<b>Target:</b>	IL12B	<b>Reactivity:</b>	Human		

**Additional Information**

<b>Conjugate:</b>	Biotinylated	<b>Uniprot:</b>	P29460
<b>Isotype:</b>	IgG1		

**Protein Information**

<b>Background:</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Synonyms:</b>	IL-12B;CLMF p40;IL-12 subunit p40;NKSF2
<b>Common Name:</b>	AK101,AK-101
<b>Recommended Dilutions:</b>	ELISA 1:5000-10000
<b>Applications:</b>	ELISA
<b>Usage:</b>	Research use only
<b>Formulation:</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage:</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.