

Anti-SELPLG (Cepeprubart Biosimilar) Antibody

HDBS0499

Description

This Anti-SELPLG (Cepeprubart Biosimilar) Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDBS0499
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Category:	Biosimilar Antibody
Synonyms:	selectin P ligand, selectin P glycoprotein ligand, PSGL-1, CD162, SELPLG
Applications:	ELISA Flow Cyt
Conjugation:	Unconjugated
Reactivity:	Human
Common Name:	cepeprubart

Antibody Data

Isotype:	Human IgG2(P193T, V282M, K322A, S378A) – kappa
Uniprot ID:	Q14242
Host Species:	Chimeric Humanized
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Target:	CD162

Preparation & Storage

Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. Store Bradford Reagent at Room Temperature for 1 Year.
Recommended Dilutions:	ELISA 1:5000-10000, Flow Cyt 1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol