

Anti-C5 (Eculizumab Biosimilar) Antibody

HDBS0533

Description

This Anti-C5 (Eculizumab Biosimilar) Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDBS0533
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Category:	Biosimilar Antibody
Synonyms:	C5D, C5a, C5b, ECLZB, CPAMD4
Applications:	N/A
Conjugation:	Unconjugated
Reactivity:	Human
Common Name:	5G1-1, 5G1.1, Soliris

Antibody Data

Isotype:	IgG2
Uniprot ID:	P01031
Host Species:	N/A
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Target:	C5

Preparation & Storage

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions: N/A

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol