

HDBS0589

Product Information

Product SKU:	HDBS0589	Host Species:	Whole mAb	Size:	100 µg
Target:	IL33	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	O95760
Isotype:	IgG4		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	DVS27;IL1F11;NF-HEV;NFEHEV;C9orf26
Common Name:	REGN-3500;REGN3500;REGN 3500;SAR-440340;SAR440340;SAR 440340
Recommended Dilutions:	-
Applications:	-
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.