

HDBS0691

Product Information

Product SKU:	HDBS0691	Host Species:	Whole mAb	Size:	100 µg
Target:	IL5	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	P05113
Isotype:	IgG4		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.		
Synonyms:	EDF;TRF;IL-5;Interleukin-5		
Common Name:	CDP-835;CDP835;CDP 835;CEP-38072;CEP38072;CEP 38072;CINQAERO;CINQAIR;Cinquil;CTx-55700;CTx55700;Ctx 55700;DCP-835;DCP835;DCP 835;SCH 55700;SCH55600;SCH-55700;TRFK-5;TRFK 5;TRFK5		
Recommended Dilutions:	-		
Applications:	-		
Usage:	Research use only		
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.		
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.		