

HDBS0758

---

**Product Information**

|                     |          |                      |           |              |        |
|---------------------|----------|----------------------|-----------|--------------|--------|
| <b>Product SKU:</b> | HDBS0758 | <b>Host Species:</b> | Whole mAb | <b>Size:</b> | 100 µg |
| <b>Target:</b>      | IL17A    | <b>Reactivity:</b>   | Human     |              |        |

---

**Additional Information**

|                   |              |                 |        |
|-------------------|--------------|-----------------|--------|
| <b>Conjugate:</b> | Unconjugated | <b>Uniprot:</b> | Q16552 |
| <b>Isotype:</b>   | IgG1         |                 |        |

---

**Protein Information**

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** IL17;CTLA8;IL-17;ILA17;CTLA-8;IL-17A

**Common Name:** TBC

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.