

Biotinylated Anti-IL17A (Bimekizumab biosimilar) mAb

HDBS0851



Product Information

| | | | | | |
|---------------------|----------|----------------------|-----------|--------------|--------|
| Product SKU: | HDBS0851 | Host Species: | Whole mAb | Size: | 100 µg |
| Target: | IL17A | Reactivity: | Human | | |

Additional Information

| | | | |
|-------------------|--------------|-----------------|--------|
| Conjugate: | Biotinylated | Uniprot: | Q16552 |
| Isotype: | IgG1 | | |

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IL17;CTLA8;IL-17;ILA17;CTLA-8;IL-17A

Common Name: Bimzelx;CDP-4940;CDP4940;UCB-4940;UCB4940

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.