

Anti-CD94 (Dibotatug Biosimilar) Biotinylated Antibody

HDBS0892

Description

This Anti-CD94 (Dibotatug Biosimilar) Biotinylated Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDBS0892
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Category:	Biosimilar Antibody
Synonyms:	KP43, KLRD1
Applications:	ELISA Flow Cyt
Conjugation:	Biotinylated
Reactivity:	Human
Common Name:	DR-01, DR01, DR 01

Antibody Data

Isotype:	IgG1
Uniprot ID:	Q13241
Host Species:	N/A
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Target:	CD94

Preparation & Storage

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Store Bradford Reagent at Room Temperature for 1 Year.

Recommended Dilutions: ELISA 1:5000-10000, Flow Cyt 15µg/ml

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol