

**Biotinylated Anti-SNCA (Exidavnemab biosimilar)
mAb**

HDBS0919



Product Information

Product SKU:	HDBS0919	Host Species:	Whole mAb	Size:	100 µg
Target:	SNCA	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P37840
Isotype:	IgG4		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	PD1;NACP;PARK1;PARK4
Common Name:	PR-1808164;PR1808164;ABBV-0805;ABBV0805;ABBV805;ABBV-805;BAN-0805;BAN0805;BAN805;BAN-805;PD-0805;PD0805;PD-805;PD805
Recommended Dilutions:	-
Applications:	-
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.