

**Biotinylated Anti-CD47 (Ligufalimab biosimilar)  
mAb**

**HDBS0970**



**Product Information**

<b>Product SKU:</b>	HDBS0970	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	CD47	<b>Reactivity:</b>	Human		

**Additional Information**

<b>Conjugate:</b>	Biotinylated	<b>Uniprot:</b>	Q08722
<b>Isotype:</b>	IgG4		

**Protein Information**

<b>Background:</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Synonyms:</b>	CD47;IAP;OA3;MER6
<b>Common Name:</b>	AK-117;AK117;AK 117
<b>Recommended Dilutions:</b>	-
<b>Applications:</b>	-
<b>Usage:</b>	Research use only
<b>Formulation:</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage:</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.