

Biotinylated Anti-Her2 (Pertuzumab biosimilar) mAb

HDBS1032



Product Information

| | | | | | |
|---------------------|----------|----------------------|-----------|--------------|--------|
| Product SKU: | HDBS1032 | Host Species: | Whole mAb | Size: | 100 µg |
| Target: | Her2 | Reactivity: | Human | | |

Additional Information

| | | | |
|-------------------|--------------|-----------------|--------|
| Conjugate: | Biotinylated | Uniprot: | P04626 |
| Isotype: | IgG1 | | |

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: NEU;NGL;HER2;TKR1;CD340;HER-2;VSCN2;MLN 19;MLN-19;c-ERB2;c-ERB-2;HER-2/neu;p185(erbB2)

Common Name: R-1273;R1273;R 1273;RG-1273;RG 1273;RG1273;RO-4368451;RO4368451;RO 4368451

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.