

Product Information

Product SKU:	HDBS1096	Host Species:	Whole mAb	Size:	100 µg
Target:	PD-L1	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	Q9NZQ7
Isotype:	IgG4		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1
Common Name:	Cejemly;CS-1001;CS1001;CS 1001;Sugelimab;WBP-3155;WBP3155;WBP 3155;Zegema;Zegemet;Zegemex;ZEJIEMEI
Recommended Dilutions:	-
Applications:	-
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.