

**Biotinylated Anti-IL17A (Vunakizumab biosimilar)
mAb**

HDBS1148



Product Information

Product SKU:	HDBS1148	Host Species:	Whole mAb	Size:	100 µg
Target:	IL17A	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	Q16552
Isotype:	IgG1		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	IL17;CTLA8;IL-17;ILA17;CTLA-8;IL-17A
Common Name:	SHR 1314;SHR-1314;SHR1314
Recommended Dilutions:	-
Applications:	-
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.