

# Biotinylated Anti-IL36RN (Spesolimab biosimilar) mAb

HDBS1179



## Product Information

<b>Product SKU:</b>	HDBS1179	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	IL36RN	<b>Reactivity:</b>	Human		

## Additional Information

<b>Conjugate:</b>	Biotinylated	<b>Uniprot:</b>	Q9UBH0
<b>Isotype:</b>	IgG1		

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** FIL1;FIL1D;IL1F5;IL1L1;PSORP;IL1HY1;IL1RP3;IL36RA;IL-36Ra;PSORS14;FIL1(Delta)

**Common Name:** BI 655130;BI655130;BI-655130;SPEVIGO

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.