

Anti-MMAE [8B4] Biotinylated Antibody

HDBS1191

Description

This Anti-MMAE [8B4] Biotinylated Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDBS1191
Contents:	10µg Bradford Reagent: 1 vial (2ml)
Category:	Payload Antibody
Synonyms:	N.A.
Clone:	Monoclonal
Clone ID:	8B4
Applications:	ELISA
Conjugation:	Biotinylated
Reactivity:	N.A.

Antibody Data

Isotype:	Rabbit IgG
Uniprot ID:	N.A.
Host Species:	Rabbit
Purification:	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Target:	MMAE

Preparation & Storage

Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Store Bradford Reagent at Room Temperature for 1 Year.
Recommended Dilutions:	ELISA 1:5000-10000
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol