

## Anti-DLL3 (tarlatamab without CD3 Biosimilar) Antibody

HDBS1198

### Description

---

This Anti-DLL3 (tarlatamab without CD3 Biosimilar) Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	HDBS1198
<b>Contents:</b>	100µg Bradford Reagent: 1 vial (2ml)
<b>Category:</b>	Biosimilar Antibody
<b>Synonyms:</b>	SCDO1
<b>Applications:</b>	<b>Flow Cyt</b>
<b>Conjugation:</b>	Unconjugated
<b>Reactivity:</b>	Human
<b>Common Name:</b>	AMG-757 without CD3

### Antibody Data

---

<b>Isotype:</b>	Human IgG1 – kappa
<b>Uniprot ID:</b>	Q9NYJ7
<b>Host Species:</b>	Humanized
<b>Formulation &amp; Reconstitution:</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Target:</b>	DLL3

## Preparation & Storage

---

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Store Bradford Reagent at Room Temperature for 1 Year.

**Recommended Dilutions:** Flow Cyt 1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol