

Anti-DLL3 and CD3 (tarlatamab Biosimilar) BsAb

HDBS1204

Description

This Anti-DLL3 and CD3 (tarlatamab Biosimilar) BsAb is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDBS1204
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Category:	Biosimilar Antibody
Synonyms:	DLL3 and CD3
Applications:	ELISA
Conjugation:	Unconjugated
Reactivity:	Human
Common Name:	AMG-757

Antibody Data

Isotype:	Human IgG1(R292C, N297G, V302C)
Uniprot ID:	Q9NYJ7 and P07766
Host Species:	Humanized
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Target:	DLL3 and CD3

Preparation & Storage

Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. Store Bradford Reagent at Room Temperature for 1 Year.
Recommended Dilutions:	ELISA 1:5000-10000
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol