

## Anti-FCGR3A (AFM13\_VH\_2 without CD30 Biosimilar) Antibody

HDBS1206

### Description

---

This Anti-FCGR3A (AFM13\_VH\_2 without CD30 Biosimilar) Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                      |                                                |
|----------------------|------------------------------------------------|
| <b>SKU:</b>          | HDBS1206                                       |
| <b>Contents:</b>     | 100µg<br>Bradford Reagent: 1 vial (2ml)        |
| <b>Category:</b>     | Biosimilar Antibody                            |
| <b>Synonyms:</b>     | FCGR3                                          |
| <b>Applications:</b> | <a href="#">ELISA</a> <a href="#">Flow Cyt</a> |
| <b>Conjugation:</b>  | Unconjugated                                   |
| <b>Reactivity:</b>   | Human                                          |
| <b>Common Name:</b>  | LSIV21                                         |

### Antibody Data

---

|                                          |                                                                                                                                                                                                  |
|------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Isotype:</b>                          | Human IgG1(L234A, L235A) – Kappa                                                                                                                                                                 |
| <b>Uniprot ID:</b>                       | P08637                                                                                                                                                                                           |
| <b>Host Species:</b>                     | Homo sapiens                                                                                                                                                                                     |
| <b>Formulation &amp; Reconstitution:</b> | Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| <b>Target:</b>                           | FCGR3A                                                                                                                                                                                           |

## Preparation & Storage

---

|                                           |                                                                                                                                                                                                                                                                                                                            |
|-------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Storage:</b>                           | <p>Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.</p> <p>Store Bradford Reagent at Room Temperature for 1 Year.</p> |
| <b>Recommended Dilutions:</b>             | ELISA 1:5000-10000, Flow Cyt 1:100                                                                                                                                                                                                                                                                                         |
| <b>Protein Quantification (Optional):</b> | <p>To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol</p>                                                              |