

Anti-HEL Human IgG1 (L234A/L235A)-Kappa Isotype control Antibody

HDBS1207

Description

This Anti-HEL Human IgG1 (L234A/L235A)-Kappa Isotype control Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDBS1207
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Category:	Biosimilar Antibody
Synonyms:	HEL
Applications:	N/A
Conjugation:	Unconjugated
Reactivity:	Human
Common Name:	Hen egg Lysozyme

Antibody Data

Isotype:	Human IgG1(L234A, L235A) – Kappa
Uniprot ID:	N/A
Host Species:	Chimeric
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Target:	HEL

Preparation & Storage

Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. Store Bradford Reagent at Room Temperature for 1 Year.
Recommended Dilutions:	N/A
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol