

Anti-LIGHT (quisovalimab Biosimilar) Antibody

HDBS1213

Description

This Anti-LIGHT (quisovalimab Biosimilar) Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	HDBS1213
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Category:	Biosimilar Antibody
Synonyms:	TNFSF14, HVEM, UNQ391/PRO726
Applications:	ELISA Flow Cyt
Conjugation:	Unconjugated
Reactivity:	Human
Common Name:	CERC-002, MDGN-002, SAR-252067, F-19, AEVI-002

Antibody Data

Isotype:	Human IgG4(L235E) – kappa
Uniprot ID:	O43557
Host Species:	Homo sapiens
Formulation & Reconstitution:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Target:	LIGHT

Preparation & Storage

Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. Store Bradford Reagent at Room Temperature for 1 Year.
Recommended Dilutions:	ELISA 1:5000-10000, Flow Cyt 1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol