

VivoGenie SARS-CoV-2 Spike RBD [2381] Low Endotoxin Antibody

IVMB0233

Description

This VivoGenie SARS-CoV-2 Spike RBD [2381] Low Endotoxin Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	IVMB0233
Contents:	100ug, 500ug Bradford Reagent: 1 vial (2ml)
Synonyms:	COV2-2381, SARS-CoV-2 Spike RBD Antibody, Receptor Binding Domain Antibody
Category:	Recombinant Monoclonal Antibody
Target:	SARS-CoV-2 RBD
Clone:	2381
Isotype:	Human IgG1
Applications:	ELISA
Specificity:	Anti-SARS-CoV-2 Spike RBD, clone 2381, specifically targets an epitope on the SARS-CoV-2 spike protein receptor-binding domain (RBD).

Antibody Data

Reactivity:	SARS-CoV-2, Virus
Host species:	-
Expression Host:	HEK-293 Cells
Immunogen:	Sequenced from Human survivors of COVID-19 (SARS-CoV-2)

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Product concentration:	≥1.0 mg/ml
Endotoxin Level:	-
Purity:	≥90% Monomer by analytical SEC and SDS-Page
Formulation:	This recombinant monoclonal antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.2 - 7.4 with no carrier protein, potassium, calcium or preservatives added. Due to inherent biochemical properties of antibodies, certain products may be prone to precipitation over time. Precipitation may be removed by aseptic centrifugation and/or filtration.

Preparation & Storage

Storage:	This antibody may be stored sterile as received at 2-8°C for up to one month. For longer term storage, aseptically aliquot in working volumes without diluting and store at ≤ -70°C. Avoid Repeated Freeze Thaw Cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Shipping:	Standard Overnight on Blue Ice.
Preparation:	Recombinant antibodies are manufactured in an animal free facility using only in vitro protein free cell culture techniques and are purified by a multi-step process including the use of protein A or G to assure extremely low levels of endotoxins, leachable protein A or aggregates.
Recommended Dilution Buffer:	-
Recommended Usage:	-
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol