

Anti-Human IL 12/23 (Briakinumab) Dylight® 488 Biosimilar

IVMB0511

Description

This Anti-Human IL 12/23 (Briakinumab) Dylight® 488 Biosimilar is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU: IVMB0511

Contents: 100ug
Bradford Reagent: 1 vial (2ml)

Synonyms: IL-12p40, Interleukin 12, Interleukin 23, IL12, IL23, IL-12, IL-23

Category: Biosimilar Recombinant Human Monoclonal Antibody

Target: IL-12/IL-23 p40

Clone: ABT-874

Isotype: Human IgG1 λ

Applications: **FA** **FC** **IF**

Specificity: This non-therapeutic biosimilar antibody uses the same variable region sequence as the therapeutic antibody Briakinumab. Briakinumab recognizes both human IL12 and IL23 via IL-12/23p40. This product is for research use only.

Antibody Data

Reactivity: Human

Host species: Human

Expression Host: HEK-293 Cells

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Immunogen: This antibody was produced by phage display technology.

Product concentration: 0.2 mg/ml

Endotoxin Level: -

Purity: -

Formulation: This DyLight® 488 conjugate is formulated in 0.01 M phosphate buffered saline (150 mM NaCl) PBS pH 7.4, 1% BSA and 0.09% sodium azide as a preservative.

Preparation & Storage

Storage: This DyLight® 488 conjugate is stable when stored at 2-8°C. Do not freeze. Store Bradford Reagent at Room Temperature for 1 Year.

Shipping: Next Day 2-8°C

Preparation: -

Recommended Dilution Buffer: -

Recommended Usage:	Application	Recommended Usage
	FC	The suggested concentration for Briakinumab biosimilar antibody for staining cells in flow cytometry is ≤ 1.0 µg per 10 ⁶ cells in a volume of 100 µl. Titration of the reagent is recommended for optimal performance for each application.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol